



## Tanascan 3D scanning system



- Scan all shapes in 3D
- Whatever material or color.
- Pause-proceed feature during scanning.
- **No post processing required** : scan and modify directly
- Image speed of 30 pps (the fastest on the market ; no tracking loss)
- Your tablet displace the result directly
- CAD compatible format

## Modification : 3 solutions supported by all scanners :

### TANASCULT : (supplied with our scan)

- All your 3D Shape scan be modified : Seat braces, mattresses, orthosis arms, etc.. unlimited possibilities.
- Simple intuitive and rapid software
- With the usual CAD tools.

### TANASCULT PRO

- The same simplicity but with added functions
- Boolean operation for removal of positive shape from soft foam

### TANACAD : Orthotists and Prosthetists have their own professional interface on SolidWork

- **Starting with measurements or 3Dfiles**, create your own positive or negative shapes.
- The simplicity of a professional interface **combined with the most powerful 3D software**.
- **Step by step guidance** during your modification
- **No shape limits**, either negative or positive
- Can be exported to a **3Dprinter**
- The main function are : 3D shape import (STL, PLY,.....), automatic searching in your database for the best adapted shapes, unlimited access to your photographs, X-rays and 3D files, design your own products (cutouts, joints and all 3D projects.)

